



High-level  
staining technology  
Now  
Made in Germany

# Leica Autostainer XL

Programmable Multiple  
Routine Stainer  
For High Flexibility and Reliability  
in Histology and Cytology

**Leica**  
MICROSYSTEMS

# Flexibility, High Specimen Throughput and Maximum Quality and Safety

The Leica Autostainer XL for routine staining in the histopathology laboratory has set new standards: The instrument having only been on the market for a few years – already more than 2,000 labs worldwide have decided on the Autostainer XL principle of innovation plus flexibility plus reliability, thus converting the instrument into one of the most successful within its class. With our new instrument generation – now made in Germany – we have increased the standard in histological routine staining even further. Very important: the Leica Autostainer XL can be used for H&E routine staining alone or simultaneously with other staining protocols. Without exception, user of the instrument confirm that the tissue sections adhere firmly to the slides and staining results are exceptionally reproducible. Because of its simple and logical operation daily work with the Autostainer XL is especially easy and efficient.

The Leica Autostainer XL's precision-made, trouble-free, stainless steel robot arm moves swiftly and quietly on the x, y and z direction motion axes to process the slide racks.



The microprocessor controlled Leica Autostainer XL has 15 programs and is designed to be simple to use and operate.



The Leica Autostainer XL has an innovative continuous load and unload drawer access and robot arm system that lets you continue slide processing without stopping the instrument's operation.



## Leica CV5030 Glass Coverslipper

**A glass coverslipper flexible enough to meet the needs of your laboratory.**

The Leica CV5030 Glass Coverslipper easily adapts to the specific requirements of varying samples such as histology sections, cytology smears, and monolayers. The CV5030 economically produces slides of superior optical quality and reliable long-term storage. The acceptance of common brands of mounting media and different sizes of coverglass allows the CV5030 to conform to your current laboratory environment. The simple operation and advanced engineering of the CV5030 will provide years of reliable use.

## Leica Autostainer XL – at one glance...

- Compact, space saving H&E stainer for routine applications
- High specimen throughput
- Simultaneous realization of various different staining protocols
- Exact incubation times
- Continuous slide unload / reload function without having to interrupt a staining cycle or open the lid
- Simple, menu-driven programming
- Incubation times and sequence of use of reagent stations freely programmable
- Integrated oven for optimal slide drying
- Reagent containers can be exchanged quickly and easily
- Integrated fume extraction system with activated charcoal filter
- Minimized user exposure to hazardous reagent fumes
- Easy-to-clean and resistant surfaces made out of polyester epoxy resin and stainless steel.






Fully automated staining  
with the  
**Leica Autostainer XL**  
and  
fully automated coverslipping  
with the  
**Leica CV5030**



**... and now completely new:**

**Autostainer XL Workstation**  
Autostainer XL Workstation




**A totally automated system  
featuring walk-away convenience**


As soon as the staining procedure is completed, an innovative robot arm transfers the slide racks from the stainer to the coverslipper. Absolutely no user intervention is required.




Chronic  
progreidant  
hepatitis;  
van Gieson



Acute cardiac  
infarction  
(24 hours old);  
H&E



Micro embolus  
in the heart;  
H&E



Lid – Covers the entire staining operations area providing user protection against reagent fumes. Gas pressure spring prevents sudden closure of the lid.

Fume extraction system – with replaceable activated charcoal filter, which absorbs hazardous fumes.

Robot arm – With “x, y and z” direction motion picks up and places racks.

Quick reference guide

At least 18 reagent stations – each one with a capacity of 450 ml. Lids prevent reagents from evaporating. Easy to lift out reagent containers.

Rack Exit Door – Allows for continuous slide unloading without interrupting the staining process.

Rack load and unload door opening indicators.

5 wash stations – with flow control.

Rack Entry Door – Allows for continous slide loading without interrupting the staining process.


Slide rack – Each rack holds 30 slides in vertical position. Can process up to 11 racks at one time. Designed to minimize carry over of reagents.

Slide drying heating station – Temperature range OFF/30° to 65° C.

Computer display – 20 characters, 4 line LCD with LED backlighting.

Computer controlled keyboard – Membrane keyboard with four function keys, four cursor keys, plus LEDs for power, processing, alarm and unload indicators.

Dimensions:  
Width: 1090 mm  
Depth: 670 mm  
Height: 510 mm  
Weight: 65 kg



# Leica Autostainer XL – Technical Specifications

Specimen slide throughput:.....at least 200 specimen slides per hour (depending on the selected program – up to 600 slides per hour)	Oven:.....1
Loading capacity:.....11 slide racks	Oven chamber temperature:.....off or 30 °C to 65 °C
Loading capacity per slide rack:.....30 specimen slides	Incubation time setting:.....from 0 sec. up to 99 min, 59 sec.
Total number of processing stations:.....26	Load / unload stations:.....1 each
Reagent stations:.....at least 18	Permanent memory capacity:.....15 programs, up to 25 program steps each
Reagent container volume:.....450 ml	Dimensions (W x D x H):.....109 cm x 67 cm x 51 cm 42.2 in. x 26 in. x 20 in.
Number of wash stations:.....max. 5	Weight:.....65 kg/143 lbs.
	Voltages:.....110 V - 120 V, 50 Hz - 60 Hz 230 V - 240 V, 50 Hz - 60 Hz

## Leica Autostainer XL Workstation

Instrument type:.....Automated Stainer (Leica Autostainer XL) and Coverslipper (Leica CV5030), modular, connected through Leica Transfer Station TS5015	Supply voltage:.....see individual instruments (TS5015 Transfer Station is fully dependent on CV5030, electrical and control)
Sample types:.....Histology and Cytology specimen	Dimensions (W x D x H):.....165 cm 67 x cm x 55 cm 65 in. x 26 in. x 22 in.
Fume handling:.....Activated carbon filter and fume extraction hose system	Dry weight:.....approx. 164 kg approx. 361 lbs.

For details on individual module characteristics please refer to the instrument data information.

Conforms to these safety and quality standards: UL/c-UL and CSA. Highly sophisticated design, manufacturing and quality assurance procedures – in compliance with DIN EN ISO 9001 – assure maximum quality and reliability.

Country of origin: Federal Republic of Germany.

Wide range of accessories on request. Technical specifications subject to change without prior notice.

## For further information please contact your local Leica sales company

<b>Australia:</b>	Gladesville/NSW	Phone +61 2 9897 9700	Fax +61 2 9817 8358
<b>Austria:</b>	Vienna	Phone +43 1 486 80 50	Fax +43 1 486 80 50 30
<b>Canada:</b>	Richmond Hill/Ontario	Phone +1 905 762 2000	Fax +1 905 762 8937
<b>Denmark:</b>	Herlev	Phone +45 4454 0101	Fax +45 4454 0111
<b>France:</b>	Rueil-Malmaison Cedex	Phone +33 1 4732 8585	Fax +33 1 4732 85 86
<b>Germany:</b>	Bensheim	Phone +49 6251 136 0	Fax +49 6251 136 155
<b>Italy:</b>	Milan	Phone +39 0257 486.1	Fax +39 0257 40 3273
<b>Japan:</b>	Tokyo	Phone +81 3 5435 9603	Fax +81 3 5435 9615
<b>Korea:</b>	Seoul	Phone +82 2 514 6543	Fax +82 2 514 6548
<b>Netherlands:</b>	Rijswijk	Phone +31 70 4132130	Fax +31 70 4132139
<b>People's Republic of China:</b>	Hong Kong	Phone +852 2 564 6699	Fax +852 2 564 4163
<b>Portugal:</b>	Lisbon	Phone +351 1 388 9112	Fax +351 1 385 4668
<b>Singapore:</b>	Singapore	Phone +65 6779 7823	Fax +65 6773 0628
<b>Spain:</b>	Barcelona	Phone +34 93 494 9530	Fax +34 93 494 9532
<b>Sweden:</b>	Sollentuna	Phone +46 8 6254 545	Fax +46 8 6254 510
<b>Switzerland:</b>	Glattbrugg	Phone +41 1 809 34 34	Fax +41 1 809 34 44
<b>United Kingdom:</b>	Milton Keynes	Phone +44 1 908 246246	Fax +44 1 908 609992
<b>USA:</b>	Bannockburn/Illinois	Phone +1 847 405 0123	Fax +1 847 405 0164



...and more than 100 national distributors.

Leica Microsystems Nussloch GmbH  
Heidelberger Strasse 17-19  
D-69226 Nussloch

Telephone: (06224) 143-0

Fax: (06224) 143 200

e-mail: [histo\\_info@leica-microsystems.com](mailto:histo_info@leica-microsystems.com)

[www.leica-microsystems.com](http://www.leica-microsystems.com)

