

Warm Stages (Click here to link to Website)

Warm stages are plates that will heat a microscope slide or petri dish up to 60°C. A controller is required to heat and control the temperature. The MC60 is suitable for applications up to 400x total magnification. A DC60 is required in high magnification applications up to 1000X. An MC60D is needed to simultaneously control two warm stages.

Please Note: Warm stages marked ** are not compatible with the DC60 controller

Standard Warm Stages

13000	TH 60-6	Warm Stage 75 x 50 x 1.4mm with 6mm aperture	£126.73
13001	TH 60-68	Warm Stage 75 x 50 x 1.4mm with 6 x 8 mm aperture	£126.73
13002	THL 60-8	Warm Stage 100 x 100 x 1.4mm with 8mm aperture	£171.36
13003	THL 60-12	Warm Stage 100 x 100 x 1.4mm with 12mm aperture	£171.36
13004	THL 60-16	Warm Stage 100 x 100 x 1.4mm with 16mm aperture	£171.36
13005	THL 60-12x25	Warm Stage for viewing microscope slides with 12 x 25mm aperture	£230.26

Nikon Specific

13006	THN 60-10	Warm Stage for NIKON Diaphot and TE series inverted microscope with 10mm aperture	£255.25
13007	THN 60-16	Warm Stage for NIKON Diaphot and TE series inverted microscope with 16mm aperture	£255.25
13008	THN 60-15x20	Warm Stage for NIKON Diaphot and TE series inverted microscope with 15 x 20mm aperture	£282.03
13009	TH 60-SMZ**	Warm Stage 205 x 155 x 2mm for NIKON Stereomicroscope with 16mm aperture	£342.72

Olympus Specific

13010	THO 60-10	Warm Stage for OLYMPUS IMT2, IX 70/50 inverted microscope with 10mm aperture	£255.25
13011	THO 60-16	Warm Stage for OLYMPUS IMT2, IX 70/50 inverted microscope with 16mm aperture	£255.25
13012	THO 60-15x20	Warm Stage for OLYMPUS IMT2, IX 70/50 inverted microscope with 15 x 20mm aperture	£282.03
13013	TH 60-SZH**	Warm Stage 205 x 155 x 2mm for OLYMPUS SZH and SZX Stereomicroscope with 16mm aperture	£342.72

Zeiss Specific

13014	THA 60-16	Warm Stage for ZEISS Axiovert inverted microscope with 16mm aperture (Axiovert 200M)	£282.03
13015	THA 60-20	Warm stage for Zeiss Axiovert inverted microscope with 20mm aperture	£282.03
13016	TH 60-SZA**	Warm Stage for ZEISS Stereozoom microscope with 16mm aperture	£340.93
13017	THA 100	Warm table with inbuilt stage for Zeiss Axiovert 100 with 16 mm aperture	£367.71
13018	THA 25	Warm table with inbuilt stage for Zeiss Axiovert 25	£399.84
13019	WS200-20	Warm Stage for Zeiss Stemi 1000 and 2000 Microscope	£340.93
13020	THZAV-60-16	Warm Stage for Zeiss Axiovert 100/135/200(to replace K-plate in universal stage)	£449.82
13034	THA 200**	Warm table with inbuilt stage for Zeiss Axiovert 200 with 16 mm aperture	£367.71
22222	THZ60-SK 12x25	Warm table for Axioskop 40 table (DV2339)	£289.17
22222	THZ60-ST 12x25	Warm table for AxioStar table (DV2340)	£289.17

Leica Specific

13021	THF 60-16**	Warm Stage for LEICA Fluovolt & Labovolt microscopes with 16mm aperture	£255.25
13022	THD 60-16	Warm Stage for LEICA DMIRB inverted microscope	£255.25
13023	TH 60-W80**	Warm Stage (205 x 155 x 2mm) for the Leica Wilde with 80mm base insert	£342.72
13024	TH 60-W120**	Warm Stage (205 x 155 x 2mm) for the Leica Wilde with 120mm base insert	£342.72
13025	WS200-HL**	Warm Stage for Leica stereo Transmitted Light Base(205x155mm with 10mm aperture) MZ12.5 Plate	£390.91
13026	WS110R-DMIL**	Warm stage for Leica DMIL (110mm dia.) (Table must be sent to Linkam for modification)	£355.21
13027	THDML60-8**	Warm Stage for Leica DMLA microscope with 8mm aperture	£251.68
13028	THDML60-12**	Warm Stage for Leica DMLA microscope with 12 aperture	£251.68
13029	THDML60-16**	Warm Stage for Leica DMLA microscope with 16mm aperture	£251.68
13030	THDML60-12x25**	Warm Stage for Leica DMLA microscope with 12x25 mm aperture	£289.17
13032	WS200-10-120**	Warm Stage for Leica MZ12.5 + Leica Base 10 445 387 - 10mm aperture fitted with 120mm ring	£390.91
13038	WS Leica Ergostage MC	Warm Stage for Leica Ergostage table for use with MC60 controller - Leica part number 11888186	£289.17
13039	WS Leica Ergostage DC	Warm Stage for Leica Ergostage table for use with DC60 controller - Leica part number 11888186	£289.17

Misc.

13033	TH60-VIS**	Warm Stage (205 x 155 x 2mm) with 16mm aperture for Vision Alpha microscope	£342.72
22222	THMX-60	Warm Stage (100 x 100mm) with 17x3mm aperture for MX table to suit Meiji MX4200L upright	£342.72
13040	WSAX200	Warm Stage to fit Prior Optiscan ES111 and ES11E50 XY tables	£449.82

Warm Stage Controllers

Please Note: Warm stages marked ** are not compatible with the DC60 controller

14023	MC60	Basic Warm stage Controller	£406.98
14024	DC60	Advanced warm stage controller for high magnification applications	£651.52
14030	MC60D	Dual channel warm stage controller for simultaneous control of two stages	£449.82